



... a symbol of top designs and products of the highest quality

CERTIFICATES OF QUALITY ASSURANCE



ÖKO-TEX




It guarantees the absence of harmful substances in our garments that are tested annually by the Trial Institute Hohenstein in Germany. This certificate is a guarantee of the quality and safety of our products. All garments are Latex-free.



ISO 9001:2001



This is an internationally recognized quality management system which certifies our streamlined policies, processes and producers. Our manufacturing efficiencies are aimed to achieve complete customer satisfaction. Our processes and procedures are inspected annually to assure compliance.



DECLARATION OF CONFORMITY



LIPOELASTIC® products meet the requirements of the Medical Device Directive 93/42/EEC and the Government Regulation No. 336/2004 as amended. The LIPOELASTIC® products are authorized to use the CE mark.

DESIGN OPTIONS

WITH ZIPPERS




Design with zippers at the sides.

WITHOUT ZIPPERS



Design without any zippers, ideal for the second stage of long-time compressive therapy.

COMFORT




Zippers with hook underlays allow for easy dressing.

VARIANT




Three position hooks offer the ability to adjust the garment after the swelling subsides.

HYGIENIC OPENING TYPES



All male and female garments with trouser-legs have a hygienic opening for increased comfort and accessibility.



All male and female garments without trouser-legs have a hygienic opening with hooks and loops for increased comfort and accessibility.

BANDAGES DESIGN



The front double underlay is now available in the newest design.



Trendy lace has a silicone wavelet on the reverse side. Therefore, the garment remains fixed in place.

CARE

LIPOELASTIC® compressive garments are machine or manually washable at the temperature of 30°C (86°F) max. Before washing, fasten the zip-fasteners and velcro fasteners. Do not squeeze and iron, dry in a shady and airy place. Drying at a high temperature reduces the elasticity of the material and reduces its curative effect.

CONTENT LIST OF THE PRODUCTS

FEMALE BANDAGES

VH	9
VH special	10
VH variant	11
VH body	12
VG	13
VG variant	14
VF	15
VF variant	16
VF body	17
VD	18
VD variant	19
VD body	20
VD special	21
VD special variant	22
VB	23
VB variant	24
VB body	25
TG	26
TG variant	27
TF	28
TF variant	29

TD	30
TD variant	31
TD special	32
TD special variant	33
TB	34
TB variant	35

FEMALE VESTS & BODY SUITS

MT smooth	38
MT formed	39
MT special	40
MH	41
MH exclusive	42
MH special	43
MHG	44
MHF	45
MHD	46
MHB	47
MGG	48
MGF	49
MGD	50
MGB	51

BREAST BANDAGES & STABILIZER BANDS

SI	54
SI special	55
SI formed	56
PI formed	57
PS formed	58
PI special	59
PS special	60
PI extra	61
PS extra	62
PS variant	63
PI super	64
PI exclusive	65
PI medical	66
PI standard	67

MALE BANDAGES

MTm	70
MTm variant	71
MTmL	72
MTmL variant	73
VHm	74

VHm variant	75
MGm	76
MGm variant	77

FACE BANDAGES

FM	80
FM formed	81
FM special	82
FM extra	83
FM whole face	84
PU	85


OTHERS

AP	88
AS	89
KP	90
KP special	91
SK	92
FD	93

STOCKINGS

LIPOPATIENT	96
LIPONURSE	97
LIPOTROMB	98

HOW TO MEASURE



BICEPS

Measure around the arm at the fullest part of the biceps.

BUST AND CHEST

Measure around the fullest part of the bust/ chest area while wearing usual undergarments.

UNDER BUST

Measure around the chest circumference directly under the breast.

WAIST

Measure around the smallest part of the waistline.

HIP/THIGH AREA

Stand with both feet together and measure around the fullest part of the hip/thigh area.



FEMALE BANDAGES

THE UNIQUE COLLECTION OF LIPOELASTIC®

The materials used by LIPOELASTIC® company guarantee the ideal compression with an emphasis on comfort. They guarantee elasticity and tensibility in all directions which assures stable compression and equal pressure for the treatment areas.

Our goal at LIPOELASTIC® is to assure the satisfaction of all patients. That's why we produce our garments in 3 colors and in sizes SSS through XXXXL. The bandages are produced by the highest quality sewing machines. The sewing technology applied ensures flat, strong and resistant seams. The entire the production includes up to 30 detailed steps.



VH

- USE:** Abdominoplasty, Liposuction of abdomen
- TYPES:** with zipper, without zipper or Comfort
- COLORS:** natural, black



Size	Measurements in cm		
	Under breasts	Waist	Hips
SS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128





VH special

USE: Abdominoplasty, Liposuction of abdomen
TYPES: with zippers, without zippers or Comfort
COLORS: natural, black



Size	Measurements in cm		
	Under breasts	Waist	Hips
SS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128

VH variant

USE: Abdominoplasty, Liposuction of abdomen
TYPE: variant - closure by three-position hooks
COLORS: natural, black



Size	Measurements in cm		
	Under breasts	Waist	Hips
SS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128





VH body

USE: Abdominoplasty, Liposuction of abdomen and back

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm		
	Under breasts	Waist	Hips
SS	61-66	58-63	84-88
S	67-72	64-69	89-94
M	73-79	70-76	95-101
L	80-88	77-85	102-110
XL	89-97	86-94	111-119
XXL	98-106	95-103	120-128

VG

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips and thighs

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm			
	Under breasts	Waist	Hips	Thigh
SS	61-66	58-63	84-88	48-50
S	67-72	64-69	89-94	51-55
M	73-79	70-76	95-101	56-60
L	80-88	77-85	102-110	61-65
XL	89-97	86-94	111-119	66-70
XXL	98-106	95-103	120-128	71-75





VG variant

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips and thighs

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm			
	Under breasts	Waist	Hips	Thigh
SS	61-66	58-63	84-88	48-50
S	67-72	64-69	89-94	51-55
M	73-79	70-76	95-101	56-60
L	80-88	77-85	102-110	61-65
XL	89-97	86-94	111-119	66-70
XXL	98-106	95-103	120-128	71-75

VF

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips and thighs

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm				
	Under breasts	Waist	Hips	Thigh	Above knee
SS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59





VF variant

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips and thighs

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm				
	Under breasts	Waist	Hips	Thigh	Above knee
SS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59

VF body

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips and thighs

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm				
	Under breasts	Waist	Hips	Thigh	Above knee
SS	61-66	58-63	84-88	48-50	32-34
S	67-72	64-69	89-94	51-55	35-39
M	73-79	70-76	95-101	56-60	40-44
L	80-88	77-85	102-110	61-65	45-49
XL	89-97	86-94	111-119	66-70	50-54
XXL	98-106	95-103	120-128	71-75	55-59





VD

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips, thighs and knees

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44

VD variant

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips, thighs and knees

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44





VD body

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs and knees

TYPES: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44

VD special

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs and knees

TYPES: with zippers, without zippers or Comfort

DOUBLE UNDERLAY INSIDE AND OUTSIDE OF THIGHS.

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44





VD

special variant

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips, thighs and knees

TYPE: variant - closure by three-position hooks

DOUBLE UNDERLAY INSIDE AND OUTSIDE OF THIGHS.

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44

VB

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips, thighs, knees and calves

Stabilization of calf implants.

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44





VB variant

USE: Abdominoplasty, Liposuction of abdomen, buttocks, hips, thighs, knees and calves

Stabilization of calf implants.

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44

VB body

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs, knees and calves

Stabilization of calf implants.

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61-66	58-63	84-88	48-50	32-34	27-29
S	67-72	64-69	89-94	51-55	35-39	30-32
M	73-79	70-76	95-101	56-60	40-44	33-35
L	80-88	77-85	102-110	61-65	45-49	36-38
XL	89-97	86-94	111-119	66-70	50-54	39-41
XXL	98-106	95-103	120-128	71-75	55-59	42-44





TG

USE: Liposuction of lower abdomen, buttocks, hips and thighs

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black




Size	Measurements in cm		
	Waist	Hips	Thigh
SS	58-63	84-88	48-50
S	64-69	89-94	51-55
M	70-76	95-101	56-60
L	77-85	102-110	61-65
XL	86-94	111-119	66-70
XXL	95-103	120-128	71-75

TG variant

USE: Liposuction of lower abdomen, buttocks, hips and thighs

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm		
	Waist	Hips	Thigh
SS	58-63	84-88	48-50
S	64-69	89-94	51-55
M	70-76	95-101	56-60
L	77-85	102-110	61-65
XL	86-94	111-119	66-70
XXL	95-103	120-128	71-75



TF

USE: Liposuction of lower abdomen, buttocks, hips and thighs

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black




Size	Measurements in cm			
	Waist	Hips	Thigh	Above knee
SS	58-63	84-88	48-50	32-34
S	64-69	89-94	51-55	35-39
M	70-76	95-101	56-60	40-44
L	77-85	102-110	61-65	45-49
XL	86-94	111-119	66-70	50-54
XXL	95-103	120-128	71-75	55-59

TF variant

USE: Liposuction of lower abdomen, buttocks, hips and thighs

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm			
	Waist	Hips	Thigh	Above knee
SS	58-63	84-88	48-50	32-34
S	64-69	89-94	51-55	35-39
M	70-76	95-101	56-60	40-44
L	77-85	102-110	61-65	45-49
XL	86-94	111-119	66-70	50-54
XXL	95-103	120-128	71-75	55-59

TD

USE: Liposuction of lower abdomen, buttocks, hips, thighs and knees

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44

TD variant

USE: Liposuction of lower abdomen, buttocks, hips, thighs and knees

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44



TD special

USE: Liposuction of lower abdomen, buttocks, hips, thighs and knees

TYPES: with zippers, without zippers or Comfort

DOUBLE UNDERLAY INSIDE AND OUTSIDE OF THIGHS.

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44


TD special variant

USE: Liposuction of lower abdomen, buttocks, hips, thighs and knees

TYPE: variant - closure by three-position hooks

DOUBLE UNDERLAY INSIDE AND OUTSIDE OF THIGHS.

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44





TB

USE: Liposuction of lower abdomen, buttocks, hips, thighs, knees and calves

Stabilization of calf implants.

TYPES: with zippers, without zippers or Comfort

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44

TB variant

USE: Liposuction of lower abdomen, hips, thighs, knees and calves

Stabilization of calf implants.

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58-63	84-88	48-50	32-34	27-29
S	64-69	89-94	51-55	35-39	30-32
M	70-76	95-101	56-60	40-44	33-35
L	77-85	102-110	61-65	45-49	36-38
XL	86-94	111-119	66-70	50-54	39-41
XXL	95-103	120-128	71-75	55-59	42-44





FEMALE VESTS & BODY SUITS



MT smooth

USE: Liposuction of arms, back and chest
TYPES: with zipper or Comfort
COLORS: natural, black



Size	Measurements in cm			
	Breast circumference	Under breasts	Waist	Arm
SS	72-79	61-66	58-63	26-27
S	80-87	67-72	64-69	28-29
M	88-95	73-79	70-76	30-31
L	96-103	80-88	77-85	32-34
XL	104-111	89-97	86-94	35-37
XXL	112-120	98-106	95-103	38-40

MT formed

USE: Liposuction of back, Mammoplasty (Augmentation, Breast reduction and Mastopexy)
TYPES: with zipper or Comfort
COLORS: natural, black



Size	Measurements in cm		
	Breast circumference	Under breasts	Waist
SS	72-79	61-66	58-63
S	80-87	67-72	64-69
M	88-95	73-79	70-76
L	96-103	80-88	77-85
XL	104-111	89-97	86-94
XXL	112-120	98-106	95-103





MT special


- USE:** Liposuction of upper abdomen and back, Mammoplasty (Augmentation, Breast reduction and Mastopexy)
- TYPE:** variant - closure by three-position hooks
Straining and stabilization by 4 vertical lamellas
- COLORS:** natural, black



Size	Measurements in cm		
	Breast circumference	Under breasts	Waist
SS	72-79	61-66	58-63
S	80-87	67-72	64-69
M	88-95	73-79	70-76
L	96-103	80-88	77-85
XL	104-111	89-97	86-94
XXL	112-120	98-106	95-103

MH

- USE:** Abdominoplasty, Liposuction of abdomen, back and upper arms
- TYPES:** with zipper or Comfort
- COLORS:** natural, black



Size	Measurements in cm				
	Breast circumference	Under breasts	Waist	Hips	Arm
SS	72-79	61-66	58-63	84-88	26-27
S	80-87	67-72	64-69	89-94	28-29
M	88-95	73-79	70-76	95-101	30-31
L	96-103	80-88	77-85	102-110	32-34
XL	104-111	89-97	86-94	111-119	35-37
XXL	112-120	98-106	95-103	120-128	38-40





MH exclusive

USE: Abdominoplasty, Liposuction of abdomen and back

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm			
	Breast circumference	Under breasts	Waist	Hips
SS	72-79	61-66	58-63	84-88
S	80-87	67-72	64-69	89-94
M	88-95	73-79	70-76	95-101
L	96-103	80-88	77-85	102-110
XL	104-111	89-97	86-94	111-119
XXL	112-120	98-106	95-103	120-128

MH special

USE: Abdominoplasty, Liposuction of abdomen and back, Mammoplasty (Augmentation, Breast reduction and Mastopexy)

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm			
	Breast circumference	Under breasts	Waist	Hips
SS	72-79	61-66	58-63	84-88
S	80-87	67-72	64-69	89-94
M	88-95	73-79	70-76	95-101
L	96-103	80-88	77-85	102-110
XL	104-111	89-97	86-94	111-119
XXL	112-120	98-106	95-103	120-128





MHG

USE: Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks and hips

TYPES: with zipper or Comfort

COLORS: natural, black





Size	Measurements in cm					
	Breast circumference	Under breasts	Waist	Hips	Thigh	Arm
SS	72-79	61-66	58-63	84-88	48-50	26-27
S	80-87	67-72	64-69	89-94	51-55	28-29
M	88-95	73-79	70-76	95-101	56-60	30-31
L	96-103	80-88	77-85	102-110	61-65	32-34
XL	104-111	89-97	86-94	111-119	66-70	35-37
XXL	112-120	98-106	95-103	120-128	71-75	38-40

MHF

USE: Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips and thighs

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm						
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Arm
SS	72-79	61-66	58-63	84-88	48-50	32-34	26-27
S	80-87	67-72	64-69	89-94	51-55	35-39	28-29
M	88-95	73-79	70-76	95-101	56-60	40-44	30-31
L	96-103	80-88	77-85	102-110	61-65	45-49	32-34
XL	104-111	89-97	86-94	111-119	66-70	50-54	35-37
XXL	112-120	98-106	95-103	120-128	71-75	55-59	38-40





MHD

USE: Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips, thighs and knees

TYPES: with zipper or Comfort

LOWER PART OF LEG: with zippers or without zippers

COLORS: natural, black



Size	Measurements in cm							
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf	Arm
SS	72-79	61-66	58-63	84-88	48-50	32-34	27-29	26-27
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32	28-29
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35	30-31
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38	32-34
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41	35-37
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44	38-40

MHB

USE: Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips, thighs, knees and calves

TYPES: with zipper or Comfort

LOWER PART OF LEG: with zippers or without zippers

COLORS: natural, black



Size	Measurements in cm							
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf	Arm
SS	72-79	61-66	58-63	84-88	48-50	32-34	27-29	26-27
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32	28-29
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35	30-31
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38	32-34
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41	35-37
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44	38-40





MGG

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks and hips, Mammoplasty (Augmentation, Breast reduction and Mastopexy)

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm				
	Breast circumference	Under breasts	Waist	Hips	Thigh
SS	72-79	61-66	58-63	84-88	48-50
S	80-87	67-72	64-69	89-94	51-55
M	88-95	73-79	70-76	95-101	56-60
L	96-103	80-88	77-85	102-110	61-65
XL	104-111	89-97	86-94	111-119	66-70
XXL	112-120	98-106	95-103	120-128	71-75

MGF

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips and thighs, Mammoplasty (Augmentation, Breast reduction and Mastopexy)

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm					
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee
SS	72-79	61-66	58-63	84-88	48-50	32-34
S	80-87	67-72	64-69	89-94	51-55	35-39
M	88-95	73-79	70-76	95-101	56-60	40-44
L	96-103	80-88	77-85	102-110	61-65	45-49
XL	104-111	89-97	86-94	111-119	66-70	50-54
XXL	112-120	98-106	95-103	120-128	71-75	55-59





MGD

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs and knees, Mammoplasty (Augmentation, Breast reduction and Mastopexy)

TYPES: with zipper or Comfort

LOWER PART OF LEG: with zippers or without zippers

COLORS: natural, black



Size	Measurements in cm						
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	72-79	61-66	58-63	84-88	48-50	32-34	27-29
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44

MGB

USE: Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs, knees and calves, Mammoplasty (Augmentation, Breast reduction and Mastopexy)

TYPES: with zipper or Comfort

LOWER PART OF LEG: with zippers or without zippers

COLORS: natural, black



Size	Measurements in cm						
	Breast circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	72-79	61-66	58-63	84-88	48-50	32-34	27-29
S	80-87	67-72	64-69	89-94	51-55	35-39	30-32
M	88-95	73-79	70-76	95-101	56-60	40-44	33-35
L	96-103	80-88	77-85	102-110	61-65	45-49	36-38
XL	104-111	89-97	86-94	111-119	66-70	50-54	39-41
XXL	112-120	98-106	95-103	120-128	71-75	55-59	42-44





**BREAST BANDAGES
& STABILIZER BANDS**



SI

USE: Implant Stabilizer Band for separate use or in combination with surgical bra as a stabilization support of implants following Augmentation

COLORS: white, black



Size	Measurements in cm
	Breast circumference
SS	72-79
S	80-87
M	88-95
L	96-103
XL	104-111
XXL	112-120

SI special

USE: Implant Stabilizer Band for separate use or in combination with a surgical bra as a stabilization support of implants following Augmentation

COLOR: white

Size	Measurements in cm
	Breast circumference
SS	72-79
S	80-87
M	88-95
L	96-103
XL	104-111
XXL	112-120





SI formed

USE: Implant Stabilizer Band for separate use or in combination with a surgical bra as a stabilization support of implants following Augmentation

COLORS: white, black

Size	Measurements in cm	
	Breast circumference	Under breasts
60	72-80	58-62
65	77-85	63-67
70	82-90	68-72
75	87-95	73-77
80	92-100	78-82
85	97-105	83-87
90	102-110	88-92
95	107-115	93-97
100	112-120	98-102



PI formed

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLORS: white, black

Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58-62	70-72	72-74	74-76	76-78	78-80
65	63-67	75-77	77-79	79-81	81-83	83-85
70	68-72	80-82	82-84	84-86	86-88	88-90
75	73-77	85-87	87-89	89-91	91-93	93-95
80	78-82	90-92	92-94	94-96	96-98	98-100
85	83-87	95-97	97-99	99-101	101-103	103-105
90	88-92	100-102	102-104	104-106	106-108	108-110
95	93-97	105-107	107-109	109-111	111-113	113-115
100	98-102	110-112	112-114	114-116	116-118	118-120



PS formed

USE: Augmentation

COLORS: white, black



Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58-62	70-72	72-74	74-76	76-78	78-80
65	63-67	75-77	77-79	79-81	81-83	83-85
70	68-72	80-82	82-84	84-86	86-88	88-90
75	73-77	85-87	87-89	89-91	91-93	93-95
80	78-82	90-92	92-94	94-96	96-98	98-100
85	83-87	95-97	97-99	99-101	101-103	103-105
90	88-92	100-102	102-104	104-106	106-108	108-110
95	93-97	105-107	107-109	109-111	111-113	113-115
100	98-102	110-112	112-114	114-116	116-118	118-120

PI special

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLORS: white, black



Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup A	Cup B	Cup C	Cup D	
60	58-62	72-74	74-76	76-78	78-80	
65	63-67	77-79	79-81	81-83	83-85	
70	68-72	82-84	84-86	86-88	88-90	
75	73-77	87-89	89-91	91-93	93-95	
80	78-82	92-94	94-96	96-98	98-100	
85	83-87	97-99	99-101	101-103	103-105	
90	88-92	102-104	104-106	106-108	108-110	
95	93-97	107-109	109-111	111-113	113-115	
100	98-102	112-114	114-116	116-118	118-120	





PS special

USE: Augmentation

COLORS: white, black



Size	Measurements in cm				
	Under breasts in cm	Cup A	Cup B	Cup C	Cup D
60	58-62	72-74	74-76	76-78	78-80
65	63-67	77-79	79-81	81-83	83-85
70	68-72	82-84	84-86	86-88	88-90
75	73-77	87-89	89-91	91-93	93-95
80	78-82	92-94	94-96	96-98	98-100
85	83-87	97-99	99-101	101-103	103-105
90	88-92	102-104	104-106	106-108	108-110
95	93-97	107-109	109-111	111-113	113-115
100	98-102	112-114	114-116	116-118	118-120

PI extra

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLOR: white



Size	Measurements in cm	
	Breast circumference	Under breasts
SS	72-79	58-64
S	80-87	65-71
M	88-95	72-78
L	96-103	79-86
XL	104-111	87-94
XXL	112-120	95-102





PS extra

USE: Augmentation

COLOR: white



Size	Measurements in cm	
	Breast circumference	Under breasts
SS	72-79	58-64
S	80-87	65-71
M	88-95	72-78
L	96-103	79-86
XL	104-111	87-94
XXL	112-120	95-102

PS variant

USE: Augmentation

COLOR: white



Size	Measurements in cm	
	Breast circumference	Under breasts
SS	72-79	58-64
S	80-87	65-71
M	88-95	72-78
L	96-103	79-86
XL	104-111	87-94
XXL	112-120	95-102






PI super

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLORS: white, black



Size	Measurements in cm				
	Under breasts in cm	Breast circumference			
		Cup A	Cup B	Cup C	Cup D
60	58-62	72-74	74-76	76-78	78-80
65	63-67	77-79	79-81	81-83	83-85
70	68-72	82-84	84-86	86-88	88-90
75	73-77	87-89	89-91	91-93	93-95
80	78-82	92-94	94-96	96-98	98-100
85	83-87	97-99	99-101	101-103	103-105
90	88-92	102-104	104-106	106-108	108-110
95	93-97	107-109	109-111	111-113	113-115
100	98-102	112-114	114-116	116-118	118-120

PI exclusive

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLOR: white




Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58-62	70-72	72-74	74-76	76-78	78-80
65	63-67	75-77	77-79	79-81	81-83	83-85
70	68-72	80-82	82-84	84-86	86-88	88-90
75	73-77	85-87	87-89	89-91	91-93	93-95
80	78-82	90-92	92-94	94-96	96-98	98-100
85	83-87	95-97	97-99	99-101	101-103	103-105
90	88-92	100-102	102-104	104-106	106-108	108-110
95	93-97	105-107	107-109	109-111	111-113	113-115
100	98-102	110-112	112-114	114-116	116-118	118-120



PI medical

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLORS: white, black



Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58-62	70-72	72-74	74-76	76-78	78-80
65	63-67	75-77	77-79	79-81	81-83	83-85
70	68-72	80-82	82-84	84-86	86-88	88-90
75	73-77	85-87	87-89	89-91	91-93	93-95
80	78-82	90-92	92-94	94-96	96-98	98-100
85	83-87	95-97	97-99	99-101	101-103	103-105
90	88-92	100-102	102-104	104-106	106-108	108-110
95	93-97	105-107	107-109	109-111	111-113	113-115
100	98-102	110-112	112-114	114-116	116-118	118-120

PI standard

USE: Mammoplasty (Augmentation, Breast reduction and Mastopexy)

COLORS: white, natural, black



Size	Measurements in cm					
	Under breasts in cm	Breast circumference				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58-62	70-72	72-74	74-76	76-78	78-80
65	63-67	75-77	77-79	79-81	81-83	83-85
70	68-72	80-82	82-84	84-86	86-88	88-90
75	73-77	85-87	87-89	89-91	91-93	93-95
80	78-82	90-92	92-94	94-96	96-98	98-100
85	83-87	95-97	97-99	99-101	101-103	103-105
90	88-92	100-102	102-104	104-106	106-108	108-110
95	93-97	105-107	107-109	109-111	111-113	113-115
100	98-102	110-112	112-114	114-116	116-118	118-120



MALE BANDAGES



MTm

USE: Liposuction and plastic surgeries of back and chest

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm	
	Chest circumference	Waist
SS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119

MTm variant

USE: Liposuction and plastic surgeries of back and chest

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm	
	Chest circumference	Waist
SS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119





MTmL

USE: Liposuction and plastic surgeries of back and chest

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm	
	Chest circumference	Waist
SS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119

MTmL variant

USE: Liposuction of back, chest and shoulders

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm	
	Chest circumference	Waist
SS	87-91	75-79
S	92-96	80-86
M	97-102	87-93
L	103-109	94-101
XL	110-117	102-110
XXL	118-125	111-119





VHm

USE: Abdominoplasty, Liposuction of abdomen and hips

TYPES: with zipper, without zipper or Comfort

COLORS: natural, black




Size	Measurements in cm		
	Under chest circumference	Waist	Hips
SS	87-91	75-79	82-85
S	92-96	80-86	86-90
M	97-102	87-93	91-97
L	103-109	94-101	98-104
XL	110-117	102-110	105-112
XXL	118-125	111-119	113-120

VHm variant

USE: Abdominoplasty, Liposuction of abdomen and hips

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm		
	Under chest circumference	Waist	Hips
SS	87-91	75-79	82-85
S	92-96	80-86	86-90
M	97-102	87-93	91-97
L	103-109	94-101	98-104
XL	110-117	102-110	105-112
XXL	118-125	111-119	113-120





MGm

USE: Abdominoplasty, Liposuction of abdomen, chest, back, buttocks and hips

TYPES: with zipper or Comfort

COLORS: natural, black



Size	Measurements in cm			
	Chest circumference	Waist	Hips	Thigh
SS	87-91	75-79	82-85	47-50
S	92-96	80-86	86-90	51-55
M	97-102	87-93	91-97	56-60
L	103-109	94-101	98-104	61-65
XL	110-117	102-110	105-112	66-70
XXL	118-125	111-119	113-120	71-75

MGm variant

USE: Abdominoplasty, Liposuction of abdomen, chest, back, buttocks and hips

TYPE: variant - closure by three-position hooks

COLORS: natural, black



Size	Measurements in cm			
	Chest circumference	Waist	Hips	Thigh
SS	87-91	75-79	82-85	47-50
S	92-96	80-86	86-90	51-55
M	97-102	87-93	91-97	56-60
L	103-109	94-101	98-104	61-65
XL	110-117	102-110	105-112	66-70
XXL	118-125	111-119	113-120	71-75





FACE BANDAGES



FM

USE: Face, neck and chin operations, chin implant surgery

COLORS: natural, black



Size	Measurements in cm
	Head measurement (chin – top of the head)
SS	50-53
S	54-58
M	59-63
L	64-68
XL	69-73
XXL	74-78

FM formed

USE: Face, neck and chin operations, chin implant surgery

COLORS: natural, black



Size	Measurements in cm
	Head measurement (chin – top of the head)
SS	50-53
S	54-58
M	59-63
L	64-68
XL	69-73
XXL	74-78





FM special

USE: Face, neck and chin operations, chin implant surgery, ear reconstruction

COLORS: natural, black



Size	Measurements in cm	
	Head measurement (chin – top of the head)	
SS	50-53	
S	54-58	
M	59-63	
L	64-68	
XL	69-73	
XXL	74-78	

FM extra

USE: Face, neck and chin operations, chin implant surgery

COLOR: white



Size	Measurements in cm	
	Head measurement (chin – top of the head)	
SS	50-53	
S	54-58	
M	59-63	
L	64-68	
XL	69-73	
XXL	74-78	





FM

whole face

USE: Face, chin implant surgery, burns

COLORS: natural, black

Size	Measurements in cm
	Head measurement (chin – top of the head)
SS	50-53
S	54-58
M	59-63
L	64-68
XL	69-73
XXL	74-78

PU

USE: Otoplasty

COLORS: natural, black



Size	Measurements in cm
	Head measurement (over ears)
SS	50-52
S	53-55
M	56-58
L	59-61
XL	62-64
XXL	65-67





OTHERS



AP

USE: Liposuction of armpit and arms

COLORS: natural, black



Size	Measurements in cm	
	Arm circumference	Under elbow circumference
SS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34

AS

USE: Liposuction of arms

COLORS: natural, black



Size	Measurements in cm	
	Arm circumference	Under elbow circumference
SS	26-27	21-22
S	28-29	23-24
M	30-31	25-26
L	32-34	27-28
XL	35-37	29-31
XXL	38-40	32-34





KP

USE: Abdominoplasty
TYPE: velcro fasteners
COLORS: white, black

Size	Measurements in cm	
	Waist	Binder height
SS	58-63	16, 23, 31, 39
S	64-69	16, 23, 31, 39
M	70-76	16, 23, 31, 39
L	77-85	16, 23, 31, 39
XL	86-94	16, 23, 31, 39
XXL	95-103	16, 23, 31, 39

KP special

USE: Abdominoplasty
TYPE: velcro fasteners, side reinforcement by flexible strap
COLOR: white

Size	Measurements in cm	
	Waist	Binder height
SS	58-63	16, 23, 31, 39
S	64-69	16, 23, 31, 39
M	70-76	16, 23, 31, 39
L	77-85	16, 23, 31, 39
XL	86-94	16, 23, 31, 39
XXL	95-103	16, 23, 31, 39





SK

USE: Moulds the figure, making it firmer and adding tone, holds shape perfectly, making the tummy flat and shaping the hips.

COLORS: white, natural, black

Size	Measurements in cm	
	Waist	Hips
SS	58-63	84-88
S	64-69	89-94
M	70-76	95-101
L	77-85	102-110
XL	86-94	111-119
XXL	95-103	120-128

FD

USE: Liposuction of knees

COLORS: natural, black

Size	Measurements in cm	
	Above knee	Calf
SS	32-34	27-29
S	35-39	30-32
M	40-44	33-35
L	45-49	36-38
XL	50-54	39-41
XXL	55-59	42-44





STOCKINGS



LIPOPATIENT

20 DEN socks
Disposable socks

Disposable socks with antibacterial finish are designed for patients who undergo an ambulatory aesthetic surgery procedure. The socks increase comfort of patients during the operation and prevent them from getting cold. They are packaged 50 pair per box.

FEATURES:

- Antibacterial and antifungal finish decrease the occurrence of common micro-organisms.
- Universal size fits all foot sizes.

COLORS: white, beige



LIPONURSE



150 DEN AD
Compression class 17-22 mmHg

Supportive calf stockings 150 DEN AD are appreciated by all nurses. Graded compression from the ankle upwards was carefully tested and proved to significantly help blood circulation in legs. It successfully reduces pain in tired legs and stops swelling origination. It is suitable to use these stockings as a preventative measure against varices and vein diseases. The elasticity of the material increases strength to pressure points with everyday wear.

FEATURES:

- Calf stockings produced from the innovative fibre Meryl Skinlife. This material has bacteriostatic properties – preventing the growth of bacteria odor and skin irritation. LIPONURSE maintains the skin's natural balance regardless of the degree of activity.
- Anatomically balanced compression.
- Contains MICRO CAPSULES to improve skin tissue microcirculation.

COLORS: white, beige, black

Size		
SSS	34 - 36	⇒ 158
SS	36 - 38	⇒ 158
S	38 - 40	158 - 164
M	41 - 42	164 - 170
L	43 - 45	170 - 176
XL	46 - 48	176 - 182



a - Unisex shoe size EU
b - High



LIPOTROMB

Anti-Thrombotic Stockings for Immobile Patients
Compression class 13-18 mmHg (1,7 – 2,4 kPa)

LIPOELASTIC® now offers a new products which increases the safety of plastic surgery interventions. LIPOTROMB thigh stockings for prophylaxis of thrombosis with medically correct compression are designed for immobile patients before, during and after surgery. This is also a superb item for mothers immediately before and after child delivery and for other immobile patients. The purpose of compression is to create sufficient pressure on the walls of the surface and the deep veins of the lower extremities to prevent blood stagnation and to direct the blood stream towards the heart. Unlike elastic dressings these tights create constant pressure independent of correct or incorrect binding.

These stockings are made by a special knitting technology assuring simultaneous lengthwise and crosswise elasticity. The crosswise elasticity is independent of lengthwise extension. Compressive effect weakens from the distal to the proximal section.

FEATURES:

- The stockings are permeable and pleasant to the skin.
- The stockings are made of materials resistant to oils and ointments.
- The stockings do not cause adverse effects.
- The rim variant is designed for mobile patients with the rim preventing the stocking from rolling down.
- The knit strips with letters "b" or "g" clearly mark measuring positions of the extremity girth.
- The knit numbers and sizes marked with letters facilitate product identification.
- The easily exposable toe tips facilitate circulation impact checks.
- The stockings can be washed at up to 90°.

COLOR: white



Inspection hole with a filled tip and a control opening under the toes.



Heel reinforcement.



LIPOTROMB stockings are lined with a non-sliding rim 5 cm wide with silicone grid and the logo of LIPOELASTIC®.



The circle on the silicone rim marks the correctly chosen tight size

SIZE AND MEASUREMENT POINTS TABLE								
Size	1*	2*	3	4	5	6	7	8
	SSS	SS	S	M	L	XL	XXL	XXXL
b	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29	29 - 31	30 - 33
g	42 - 47	44 - 51	48 - 55	52 - 59	55 - 67	59 - 67	63 - 70	67 - 73

* Note: Sizes 1 and 2 are not mass produced and can only be made to order!

b - ankle girth is measured in the narrowest place
g - thigh girth is measured in the place of the intended rim edge

SIZE CHOICE AND CHECK

Thanks to the sophisticated LIPOELASTIC® system you will determine the right size of the LIPOTROMB stocking in a minute. The size is knit in each piece of the tight to help to check whether the correct size has been chosen for each individual patient.

The specially developed smart rings on the silicone rim will always reassure you that you have chosen the right stocking.



On the inside



On the outside



The circle on the silicone rim

lipoelastic®

marks the correctly chosen tight size

lipoelastic®

HORIZONTAL OVAL ON THE RIM
 Incorrect size but the patient will not be harmed.

lipoelastic®

VERTICAL OVAL ON THE RIM
ATTENTION – incorrect size!
 May cause serious health problems.

LIPOELASTIC a. s.

Vrbenská 887
 757 01 Valašské Meziříčí
 CZECH REPUBLIC
 Tel.: +420 571 633 512
 Fax: +420 571 633 268
 E-mail: info@lipoelastic.com
 www.lipoelastic.com



EUROPEAN UNION
 EUROPEAN REGIONAL DEVELOPMENT FUND
 INVESTMENT IN YOUR FUTURE